Instructions for Use

380094-E

English

OC-Auto 3 Latex Reagent
OC-Auto 3 Buffer

REF V-PH18, V-PH33
REF V-PH46

INTENDED USE

For quantitative measurement of haemoglobin in faeces using automated immunochemical analysers [OC-SENSOR Micro, OC-SENSOR io].

INTRODUCTION

The amount of haemoglobin in faeces increases with illnesses that are accompanied by haemorrhagic lesions in the digestive tract, particularly in the lower digestive tract. Therefore measuring the amount of haemoglobin in faeces is an effective means of screening for early detection and treatment of colon cancer and other lower digestive tract illnesses that are accompanied by hemorrhaging⁽¹⁻³⁾. Currently, the methods used for detection of haemoglobin include immunoassay latex agglutination⁹⁾, RPHA, EIA, and immunonephelometry. OC-Auto 3 Latex Reagent is an immunoassay test reagent and used together with an automatic analyser for measurement of haemoglobin in faeces. OC-Auto 3 Latex Reagent contains the latex particles sensitized to anti-human haemoglobin A₀ (HbA₀) polyclonal antibodies for optical measurement of the latex agglutination reaction.

PRINCIPLE OF THE METHOD

The test method is based on a latex agglutination reaction. A latex reagent is prepared by sensitizing anti-human HbA₀ antibodies to polystyrene latex particles. When this reagent is mixed with the sample, the anti-human HbA₀ antibodies, which were sensitized to latex, react with the haemoglobin in the sample, and the latex aggregate is formed in the latex agglutination reaction. The change in absorbance per unit time resulting from the latex agglutination reaction is proportional to the concentration of haemoglobin in the sample. A dose response curve of the absorbance unit (OD) vs. concentration is generated using the results obtained from the calibrators. The concentration of haemoglobin in the patient sample is determined from this curve.

REAGENTS

Reagents are stable until the date printed on the label assuming the container remains unopened at a storage temperature of 2-10°C.

OC-AUTO 3 Latex Reagent7mL, 2 vials

(1 vial contains 1.4mL of latex particles suspension coated with anti-human HbA₀ rabbit IgG)

OC-AUTO 3 Buffer 200mL, 1 vial

(1 vial contains 2.38g of N-2-Hydroxyethylpiperazine-N'-2-ethanesulfonic acid (HEPES))

Note: The above reagents are sold separately.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Before performing measurement, be sure to read the instruction manual of the analyser.
- Fully prepare and adjust the analyser before measurement.
- Create a calibration curve for each measurement. Also be sure to create a new calibration curve when OC-AUTO 3 Latex Reagent from a different vial or lot is used.
- Use the Calibrator (REF V-PH51, V-PH52, V-IX50) for the creation of calibration curve.
- If the measurement results exceed the measurement range, use Diluent of Calibrator to dilute the sample and perform measurement again.
- Store all reagents at 2-10°C and avoid freezing.
- Do not use reagents that have passed their expiration date.
- The test sample may contain microorganisms. Therefore, use caution when handling. After use, all samples and other materials must be considered as medical waste and disposed of accordingly.
Example of treatment: Soak for 1 hour or longer in a sodium hypochlorite solution (available chlorine concentration 1000 ppm or greater). (Neutralize any substances that contain acids before soaking.) Alternatively, treat in an autoclave at 121°C for 20 minutes. (Do not treat any items which sodium hypochlorite have adhered to in this way.)
- Dispose of used reagents and containers as medical waste in accordance with local regulations.
- If the product is used in any way other than that specified here, the reliability of measurement results cannot be guaranteed. Be sure to follow the procedure
- A clinical diagnosis based on the measurement results must be a comprehensive judgment made by the attending physician, including factors such as clinical symptoms and other test results.

SAMPLE COLLECTION DEVICE

- Use the following designated sampling container for collecting specimens.
OC-Auto Sampling Bottle 3 (REF V-PZ25, V-PZ26)
- Collect faecal samples by scraping the surface of the stools in different areas. Collect the amount necessary to cover the groove of the probe.
- Check that the test sample has become fully suspended in the buffer solution inside the sampling bottle.
- Performance testing with the sample collection device demonstrated that samples stored at 2-10°C for 28 days have 95±14.7%, at 25°C for 7 days have 96±20.4%, for 14 days have 93±23.5%, and at 30°C for 7 days have 89±20.5%, for 14 days have 84±23.6% of haemoglobin recovery (in-house data, recovery rate shown as mean±2SD). However, the haemoglobin in some samples may undergo rapid denaturation or degradation, resulting in false negatives. Therefore, samples should be stored at 2-10°C and analysed as soon as possible.

PREPARATION OF REAGENTS

- OC-AUTO 3 Latex Reagent
Ready to use.
Before measurement, bring reagent to room temperature and invert the OC-Auto 3 Latex Reagent gently several times to assure uniform suspension. Do not mix reagent with another bottle.
- OC-AUTO 3 Buffer
Ready to use.
Before use, stabilize the buffer to room temperature(20-25°C). Do not mix buffer with another lot.
Do not mix buffer from different bottles.
Do not use precipitated buffer, precipitate of buffer may cause problems.

REQUIRED MATERIALS, NOT PROVIDED BY THE MANUFACTURER (Prepare these materials before analysis)

- Wash solution: Sodium hypochlorite 0.15% (0.10%–0.30% is acceptable)
- Purified water (for washing): Distilled or de-ionized water (1.0-10.0MQ·cm is acceptable)
- Printer Paper: Use thermal printer paper which fits the analyser.

TEST PROCEDURE

Perform measurements according to the instruction manual of the analyser.

- Set the OC-AUTO 3 Latex Reagent and OC-AUTO 3 Buffer on the designated locations in the analyser.
- Input the parameters into the analyser.
- Check water, wash solution and drain tank volumes.
- Load control materials and/or sample collection devices.
- Start analysis.
- After measurement is completed, ensure the results of controls are within the acceptable range.

RESULTS

The reaction of each sample is compared to the calibration curve that was created previously. The concentration of HbA₀ (ng/mL) is found and a positive or negative judgment is made.

LIMITATIONS

- The test has not been validated for testing of patients with haemoglobinopathies.

ANALYTICAL SPECIFICITY

Almost no reaction occurs with non-human animal haemoglobin. Almost no effect on the measurement value was found from the following co-existing substances at the indicated concentrations.

| Substance | Concentration | Substance | Concentration |
|--------------------|---------------|-------------------------------|---------------|
| bilirubin | 25mg/dL | Protein(bovine serum albumin) | 2.5g/dL |
| lipids(intralipid) | 0.6% | glucose | 4.0g/dL |
| ascorbic acid | 40mg/dL | barium sulfate | 25mg/dL |

INTERNAL QUALITY CONTROL

Each laboratory should establish a quality control program to monitor the performance of OC-AUTO 3 Latex Reagent. It is recommended to use the following materials for the quality control in your laboratory.

OC-Control LV1 (REF V-PH53) OC-Control Low (REF V-PH57)
 OC-Control LV2 (REF V-PH54) OC-Control High (REF V-PH58)
 OC-Control LV3 (REF V-PH59)

PERFORMANCE CHARACTERISTICS

- Sensitivity
When standard solution of 50ng/mL HbA₀ was measured in comparison to the blind reagent test, the change in absorbance that was obtained was 0.002 or more.
- Accuracy
When control samples of known concentration were measured, the value obtained was within ±15% of the indicated value.
- Within-run reproducibility
When the same sample was measured 10 times within the same run, the coefficient of variation (CV) for the values obtained was 10% or less.
- Measurement range
50 – 1000ng/mL (10µg/g – 200µg/g (faeces))
- Cut-off value should be set on the demand of each laboratory¹⁰⁾.

PRODUCT CODE, PRODUCT NAME & STORAGE

| Product code | Product name | Contents | Storage | Analyser |
|--------------|---|-------------|---------|-----------|
| V-PH18 | OC-AUTO 3 Latex Reagent | 2×7mL | 2-10°C | io |
| V-PH33 | OC-AUTO 3 Latex Reagent | 2×7mL | 2-10°C | Micro |
| V-PH46 | OC-AUTO 3 Buffer | 1×200mL | 2-10°C | io, Micro |
| V-PH51 | OC-Calibrator 1 | 1×3mL | 2-8°C | Micro |
| V-PH52 | OC-Calibrator 2 | 1×3mL | 2-8°C | io |
| V-IX50 | OC-Standard | 10×1mL | 2-10°C | Micro |
| V-PH53 | OC-Control LV1 | 2×5mL | 2-8°C | io, Micro |
| V-PH54 | OC-Control LV2 | 2×5mL | 2-8°C | io, Micro |
| V-PH59 | OC-Control LV3 | 2×5mL | 2-8°C | io, Micro |
| V-PH57 | OC-Control Low | 5×1mL | 2-10°C | Micro |
| V-PH58 | OC-Control High | 5×1mL | 2-10°C | Micro |
| V-PZ25 | OC-Auto Sampling Bottle 3 | 100 bottles | 1-30°C | io, Micro |
| V-PZ26 | OC-Auto Sampling Bottle 3 without barcode | 100 bottles | 1-30°C | io, Micro |

REFERENCE

- J.S.Mandel et al.: New England Journal of Medicine, 328: 1365-1371, 1993.
- S.J.Winawer et al.: Journal of National Cancer Institute, 85 (16): 1311-1318, 1993.
- N.Hiwatachi et al.: Japan Journal of Cancer Research, 84: 1110-1112, 1993.
- H.Saito et al.: International Journal of Cancer, 61: 465-469, 1995.
- B. Levin et al.: Gastroenterology, 134(5): 1570-1595, 2008.
- S. Ciatto et al.: British Journal of Cancer, 96: 218-221, 2007.
- G. Grazzini et al.: European Journal of Cancer Prevention, 13: 16-26, 2006.
- E. Fenocchi et al.: European Journal of Cancer Prevention, 15: 384-390, 2006.
- T. Takeshita et al.: Journal of Coloproctology, 38: 780-783, 1985.
- LGM Rossum et al.: British Journal of Cancer, 101: 1274-1281, 2009.

Advena Ltd.

Tower Business Centre, 2nd Flr., Tower Street, Swatar, BKR 4013 Malta

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EIKEN CHEMICAL CO., LTD.